

Nitto

Innovation for Customers

XTRATA® Perme-Roll AIR™

Thin Type Film Dressing (Handy Roll Type)

More conformable with even thinner film*

*Approx. one-fourth of XTRATA® Perme-Roll



With STRATAGEL® Technology



Less damaging to the stratum corneum
Because of our revolutionary acrylic adhesive
with liquid gel component

Better conformability with superior fit to skin surface

Extra Thin Film

The extra thin film* that stretches with less stress provides pliability. It is conformable to body contour, minimizing the strain and discomfort during application. It also reduces shear stress against clothing, which often causes accidental removal of the dressing.

*Thickness is approx. 1/4 of XTRATA® Perme-Roll



STRATAGEL® Soft Gel Adhesive

STRATAGEL® reduces irritation that is caused during application and upon removal. The soft gel adhesive fits so well to skin surface that it realizes gentle removal of the dressing with minimum damage to stratum corneum.


Conventional adhesive

Skin with damaged stratum corneum




STRATAGEL®

Stratum corneum is less damaged.



Instruction for Use



- 1 Peel the paper liner in the middle.*
 - 2 Position the adhesive surface to cover absorbent pad and remove remaining paper liners on both edges.
 - 3 Remove the transparent cover film by lifting the tabs that run in the middle.
- *For L34R05, remove the larger paper liner before positioning and the other after application.

Assortment

Product Code	Size	Roll/Box
L34R05	5cm × 10m	1
L34R10	10cm × 10m	1
L34R15	15cm × 10m	1

Intended Use

- Waterproof purpose to protect a wound area that is covered with wound dressing.
- Waterproof purpose to cover external devices such as stoma, CAPD
- Fixation of a catheter and a drain

Precautions

- Do not apply the adhesive section directly onto the wound.
- Make sure the film seals surrounding area of the wound dressing before having a shower.
- The skin should be dry prior to application to ensure adequate adhesion. Shave the application area if necessary for better fitting.
- Do not stretch the film when applying on the skin.
- Make sure to remove the transparent cover film.
- Remove the film immediately if skin damage or other symptoms occur.
- Do not reuse.
- Avoid direct sunlight and store dry and cool place.

EXTRATA® range is manufactured by Nitto Denko Corporation

Nitoms, Inc.

ANT Marketing Dept., Head Office: 7F Shinagawa Seaside Park Tower, Higashishinagawa, Shinagawa-ku, Tokyo, 140-0002, Japan

TEL: +81-3-6631-1663

<http://www.nitoms.com>